Improve business efficiency • Improve employee workplace safety • Improve patient care & management



#153093 W/Rouge - #153040

Economy Folding Screen Frame #153093

The 3 sections of the screen are attached to each other by a riveted hinge which allows the frame to smoothly fold and unfold as needed. Provides instant isolation for routine, temporary and special circumstances. Anodized aluminum frame and "accordion" hinge allows frame to fold for convenient storage. Screen panels are secured to the frame by top and bottom snap in rods.

Features:

- Available in 2, 3 or 4 section models
- Three sections measures 68"H x 40.5"W
- Requires one vinyl or fabric panel per section
- Anodized aluminum construction provides a lightweight and durable solution eliminating the possibility of paint chipping or rusting
- 5 year warranty on Frames
- 90 Days on Screen Panels



