



## **Clear PETG Bedpan Rack #303100**

Omnimed bedpan racks are constructed of lightweight acrylic that is both durable and easy to clean. Racks have pre-drilled holes that allow for convenient door or wall mounting. Units measure 11 5/8 H x 12 W x 4 3/4 D.

Construction: PETG Plastic

Size: 11<sup>5</sup>/<sub>8</sub>"H x 12"W x 4<sup>3</sup>/<sub>4</sub>"D

- Tougher than acrylic
- Will provide many years of service
- Recyclable and ECO friendly
- Easy to clean and mount