

# OMNIMED Glove Box Holder

Improve business efficiency • Improve employee workplace safety • Improve patient care & management



## Glove & Masks Station -

Omnimed glove and mask dispensers are constructed of a special 3/16" thick impact resistant acrylic polymer for increased durability. These side loading units hold one box of gloves and one box of masks. They have pre-formed keyholes for wall or door mounting but may also be used as table-top dispensers. Units measure 10"H x 10 1/2"W x 4 5/8"D. These units include a 1 year warranty.

- Constructed of 3/16" clear acrylic
- Pre-drilled for mounting to walls (hardware included)
- Measures 10"H x 10 1/2"W x 4 5/8"D
- HDDDWH

**OMNIMED**

800.257.2326 • [www.omnimedstore.com](http://www.omnimedstore.com)