

# **OMNIMED** Acrylic Suture Rack

*Improve business efficiency • Improve employee workplace safety • Improve patient care & management*

---



## **Acrylic Suture Rack - #304100**

Made of special impact resistant acrylic that's more durable than traditional acrylic. Rack is top loading and has front rail to keep boxes in place when retrieving suture package. These units include a 1 year warranty.

- Pre-drilled for mounting to walls (hardware included)
- Measures 6"H x 15½"W x 6"D
- 1 Year Warranty Included

**OMNIMED**

800.257.2326 • [www.omnimedstore.com](http://www.omnimedstore.com)