

Improve business efficiency • Improve employee workplace safety • Improve patient care & management



Deluxe Wire Glove Box Holder - #305380

Omnimed wire glove box holders are constructed of heavy gauge steel wire that is powder coated for a durable and easy to clean finish. They are ideal for use in medical, lab and food prep functions or any area where cross contamination is an issue or protective gloves are needed. These glove box holders are simple to load and accommodate up to three different glove sizes to meet staff needs. Each unit can be mounted vertically or horizontally on a variety of surfaces. These units include a limited 5 year warranty.

- Durable white powder coating on steel wire
- Easy to clean, load and mount
- Mounted vertically or horizontally
- Measures 17"H by 11½"W by 4½"D