

Improve business efficiency • Improve employee workplace safety • Improve patient care & management

**Wire Basket- #350005**

Square Wire Baskets are constructed of heavy duty powder coated steel wire and designed to be used with Omnimed's Computer Monitor Stand or Infusion Stands. These baskets are easily accessible and provide an additional storage area for supplies or equipment. They can be mounted to the stand column with a standard bracket or with the Omnimed Clamping System when multiple accessories are being attached. These units include a 1 year warranty.

- Durable powder coated steel wire
- Additional storage space for charts, supplies or equipment
- Measures 6"H by 6 ¼"W by 4 ¼"D