



## **Adjustable Kidney Stand #350315**

The Omnimed Kidney Stand is a hospital-grade, ergonomically-designed, and height-adjustable cart that delivers point-of-care technology to improve workflow and increase accuracy. Omnimed proprietary lift technology offers a 10" ergonomic height adjustment range, rugged construction, and space-saving accessibility. The Omnimed Kidney Stand is perfect in any environment – patient rooms, emergency and operating rooms, radiology labs, nurses' stations – or wherever access to patient data and images is needed. Omnimed products are proudly made in the U.S.A. Unit Includes 5 year warranty.

- Faux cherry wood surface made from easy to clean thermofoil
- Painted steel sub surface adds durability and platform to capture the modular components
- 19.75" W x 6.5" D x 2" H cubbie compartment
- Work surface dimensions 32" W x 20" D
- Height adjustment 31" to 41"
- 22" diameter star base
- Two cup holder openings



**#350320 Med Drawer**  
16" W X 9.5"D X 7" H



**#350328 Glove Holder**  
10.25"W X 3.9"D X 5.25" H



**#350331 10' Hospital Grade Cord** 6.5" D X 3" H



**#741323 Tablet Holder**



**#350325 Hand Sanitizer Holder** 3"W X 3"D X 11"H



**#350716 Extended Life Battery**



**#350327 Pencil Cup Holder**  
3" D X 4" H



**#350002 Transport Basket**  
10"H X 13"W X 4"D



**#350005 Wire Basket**  
6.25"H X 4.25"W X 4.5"D



**#350326 Sanitizer Refill**  
3" W X 3"D X 5"H