

OMNIMED Heavy Weight™ Infusion Stand

Improve business efficiency • Improve employee workplace safety • Improve patient care & management



Heavy Weight™ Infusion Stand
#741300

Heavy Weight™ Infusion Stand

#741300

Features:

- Stainless steel pole
- 4 V-hooks
- Holds up to four 5,000 cc containers
- Height adjustment up to 9 feet
- Unique “Air-Glide System” controls descent of pole when unlocked
- Heavy-duty cast aluminum 5-leg base
- 27” diameter base
- Unit Includes 3 Limited 3 year warranty

OMNIMED

800.257.2326 • www.omnimedstore.com