

# Omni Tablet Holder



## Omni Tablet Holder - #741324 - With Omni Clamp

The Omni Tablet Holder can accommodate all tablets measuring from 6" to 12 1/4" in width, 6" to 12 1/4" in height and 3/4" or less in depth. The unique tension hinged design allows the operator to easily change the angle of view and the hinges ensure that the angle does not change. No tools or knobs are required to adjust the angle. The holder also includes an Omni Clamp and bracket to mount to infusion stands, transport carts and bedrails. The clamp can fasten to square or round surfaces with a depth or diameter of 2" or less. It is made from heavy gauge steel and will provide years of durable service (5-year full warranty).



OmniMed Products...  
Made in the U.S.A.

MANUFACTURERS OF **IBeam**® PRODUCTS  
**OMNIMED** INC.®

800 Glen Avenue • Moorestown, N.J. 08057

800.257.2326 • FAX: 856.359.2249

[www.omnimedbeam.com](http://www.omnimedbeam.com)



# Omni Tablet Holder

